ZWISS CLEAN MASK BASIC

ZWISS CLEAN®
Protected by
livinguard

ANTIBACTERIAL & ANTIVIRAL PROTECT*

99.999% of viruses and bacteria are eliminated in short time

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OPTIMAL FIT

Modern design and cut according to the recommendation for personal safety equipment, with variable ear flap

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WASHABLE & RECYCABLE

100% recyclable, to be disposed of sorted by materials; washable up to 10 times at 30 degrees

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SUSTAINABLE

Sustainable production and working conditions

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LIGHTFLOW AIR

Breathable material, portable up to 12 h, even sleeping is possible

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antibacterial & antiviral textile

antibakterielles & antivirales Textil

breathable material, can be worn up to 12 h, even sleeping in it is possible

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100 % Polyamide (PA)
100 % Polyamid (PA)

*The zwissCLEAN® protected by livinguard treated material has a very high antimicrobial activity (>99.9%) against a wide range of microorganisms, including Escherichia coli as a test organism. The antimicrobial activity was tested according to AATCC 100, DIN EN ISO 20743 and DIN EN ISO 18184. The fibres contain a chemical formulation that has an antimicrobial effect through a physical reaction. In the fibre polymer this counteracts the colonisation, replication and odour development of microorganisms on the treated products.

100 % to be disposed of sorted by materials

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Instructions for the proper fitting of the zwissMASK

1. Wash and disinfect hands thoroughly (20-30 seconds with soap) before and after putting on the mask.
2. Unfold the zwissMASK.
3. Place the zwissMASK on your face and put the ear loops around your ears.
4. Check the fit on chin, mouth and nose and adjust the position of the zwissMASK and the size of the ear flap when it is necessary. The ear loops can easily be resized with scissors. Use the already pre-punched markings for this purpose.
5. The masks should be sealed airtight after removal or washed immediately.

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The oronasal mask is neither a medical product according to DIN EN 14683 nor a safety product according to DIN EN 149. The mask is not suitable for use in the health and care sector or as protective equipment.

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variable ear flap, easy to adjust the size because of the already pre-punched markings

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